

## Air overlay

Alternating air overlay pressure care mattress with quiet pump.



### Features

- 17,5" air cells
- CPR valve
- Static cells in head end
- Static mode
- Very quiet pump
- Welded seams for infection control

### Specifications

|                  |                 |
|------------------|-----------------|
| Prevention       | High Risk       |
| Treatment        | Category 2*     |
| Type             | Overlay         |
| Max. user weight | 130kg / 21st    |
| Size             | 88 x 200 x 15cm |

\*With clinical assessment  
 \*Support surface assessment requires holistic assessment and cannot be decided on risk assessment score alone (EPUAP 2015).

[VIEW MORE](#)

## Air overlay

Alternating air overlay pressure care mattress with quiet pump.

The Air Overlay mattress is designed to be used in conjunction with an underlay between the bed frame and the underside of the overlay mattress.

### Alternating pressure therapy

With seventeen 5" cells, patients are provided with a soothing alternating dynamic pressure system, which when combined with an underlay mattress, offers alternating pressure therapy for medium to high risk patients.



### Mattress underlay

Use an underlay mattress as a protective layer between the bed and overlay air mattress.



### OSKA Selector

With no clear algorithm or selection guide it can be difficult to know what mattress you should be selecting for your patients. The OSKA Selector helps carers identify if the patient is on the correct mattress.



### Infection prevention & control

OSKA's ToughWeld welded seams are melted or welded together using high frequency vibrations meaning there are no holes or gaps for ingress of body fluid or contamination.

