

# Series3 Bariatric

Active therapy support surface.



**AVAILABLE FOR RENTAL**

## CPR Valve

The flush fitting CPR push button reduces the risk of snagging and accidental mattress deflation.

## Welded seams

Eliminate stitching and support infection prevention.

## Multi-stretch, moisture vapour permeable cover

Assists pressure redistribution and microclimate management.

## Cable tidy

Elevates the cable off the floor to reduce the risk of cable damage, trips and falls.

## Cell-on-cell design

Reduces the risk of bottoming out.

## Quick connectors

Quick connectors enable creation of bespoke zero pressure zones to offload pressure from vulnerable areas.

They also allow rapid removal and replacement of damaged air cells in less than a minute with no effect on therapy.

## Digital pump\*

Simple to set up and use, with manually selected weight settings.

\*N.B. Series3 Bariatric does not have an autoweigh pump.



## Specification

Prevention	Very High Risk*
Treatment	All Categories**
Max. user weight	200kg / 31.5st
Bariatric max. user weight	400kg / 63st
Medical Device Classification	Class IIa

\*With clinical assessment. Support surface selection requires holistic patient assessment and cannot be decided on risk assessment score alone (EPUAP 2015).

\*\*Both SDTI and Unstageable need to be treated as category 3-4 until depth of damage is established.

## Mattress size

Series3 Bariatric (120 x 203 x 26cm)

## Order Code

PR13-761