

Selector

Indication	Selector*	Mattress
Complex (E.g Palliative Care)		Series5 Lateral tilt mattress
		Series6 OneSurface™ Advanced hybrid mattress
Category 3 & 4		Series3 Auto Alternating air mattress
Category 1 & 2		Series4 Non-powered hybrid mattress
Very High risk		Air Overlay Alternating air overlay mattress
		Series2 High specification foam mattress
High /Very high risk		Series I Advance High specification castellated foam pressure care mattress.
High risk		Series I High risk mattress

*With clinical assessment. Support surface selection requires holistic patient assessment and cannot be decided on risk assessment score alone (EPUAP 2015).