

Integra Shell Chair

Easy-to-manoeuvre chairs can be adapted to meet various needs of residents in nursing home environments.



Features

Ideal transfer aid, provides the user with support & comfort

A gas assisted tilt in space chair that rises and tilts users to their feet

Easy adjustment due to the ergonomically designed push handle with hand lever

Velcro attached ergonomically designed headrest, which can be easily adjusted

Lap belt is fitted as standard for additional safety

Directional wheel enables the chair to be transferred in a straight line

Good hoist access for easy moving and handling

Specifications

Max. user weight	120kg / 19st
Seat Size (W x D)	45 x 44cm
Overall Size (W x D)	69 x 89cm
Warranty	2 years

Integra Shell Chair

Easy-to-manoevre chairs can be adapted to meet various needs of residents in nursing home environments.

Used as an aid for moving & handling residents in nursing home environments; the Integra Shell Chair offers a forward tilt & tilt in space movement to aid with standing and provide the user with seated support & comfort.

Ideal transfer aid

Easy to manoeuvre high quality castors (12.5cm) which integrally rotate 360 degrees. The rear 2 castors include brakes.



Easy adjustment

Ergonomically designed push handle with hand lever for easy adjustment.



Tilt in space

A gas assisted tilt function, tilts the user forward aiding with standing up from the chair as well of offering users who are seated for long periods support & comfort.



Adjustable headrest

Velcro attached ergonomically designed headrest, which can be easily adjusted.

