




## Series1



Modular foam pressure care mattress.



### Features

Modular foam

4 way stretch covers

Ideal for use on profiling beds

Available in divan bed size

Long term cost savings to the organisation

Turn periodically to increase its life span

### Specifications

| Prevention                                | High Risk*                  |
|---|-----------------------------|
| Size                                      | 88 / 105 / 120 x 200 x 14cm |
| Max. user weight                          | 213kg / 33.5st              |
| Max. user weight<br>(120cm wide mattress) | 250kg / 40st                |

\*With clinical assessment. Support surface selection requires holistic patient assessment & cannot be decided on risk assessment score alone (EPUAP 2015).



## Series1

**Modular foam mattress provides an excellent 'High Risk' rated mattress for all patients.**

Modular foam mattress with 4 way stretch cover redistributes pressure. Turned periodically to increase its lifespan, this high risk mattress provides long term cost savings to the organisation.

### Modular foam

The modular foam cubes react independently to conform to the patient/resident's body, redistributing the interface pressure. The mattress can be turned periodically to increase its life span.



### Withstand rigorous cleaning

Where blood or body fluids are present apply up to 1% sodium hypochlorite solution (10,000ppm chlorine).



### OSKA Selector

With no clear algorithm or selection guide it can be difficult to know what mattress you should be selecting for your patients. The OSKA Selector helps carers identify if the patient is on the correct mattress.



### Pressure relieving stretch covers

The 4-way stretch covers are designed to stretch in all directions so patients & residents get the full benefit of the pressure relieving foam mattress.

