




Series1

Castellated foam pressure care mattress.



Features

- Castellated foam
- 4 way stretch covers
- Ideal for use on profiling beds
- Available in divan bed size
- Long term cost savings to the organisation
- Turn periodically to increase its life span

Specifications

Prevention	High Risk*
Max. user weight	213kg / 33.5st
Max. user weight (120cm wide mattress)	250kg / 40st

Mattress size

88 x 198 x 15cm	PR13-220
105 x 198 x 15cm	PR13-162
120 x 198 x 15cm	PR13-156

Order Code

*With clinical assessment. Support surface selection requires holistic patient assessment & cannot be decided on risk assessment score alone (EPUAP 2015).




Series1

Castellated foam mattress provides an excellent 'High Risk' rated mattress for all patients.

Castellated foam mattress with 4 way stretch cover redistributes pressure. Turned periodically to increase its lifespan, this high risk mattress provides long term cost savings to the organisation.

Castellated foam

The castellated foam cubes react independently to conform to the patient/resident's body, redistributing the interface pressure. The mattress can be turned periodically to increase its life span.



Withstand rigorous cleaning

Where blood or body fluids are present apply up to 1% sodium hypochlorite solution (10,000ppm chlorine).



OSKA Selector

With no clear algorithm or selection guide it can be difficult to know what mattress you should be selecting for your patients. The OSKA Selector helps carers identify if the patient is on the correct mattress.



Pressure relieving stretch covers

The 4-way stretch covers are designed to stretch in all directions so patients & residents get the full benefit of the pressure relieving foam mattress.

